

MECH345
MODERN INSTRUMENTATION AND
EXPENRIMENTATION


LAB REPORT #1

SARAVANA KUMAR PHILLIP
NAGARAJAN
W1409905

TASK: 1

A) Oscilloscope connected with Function Generator:

Purpose : To measure the voltage and frequency with the oscilloscope for different settings in function generator.





Frequency (Function Generator)	Amplitude (Function Generator)	Voltage (Oscilloscope) (volt/div)
500 Hz	1 (Vpp)	2.2 V
1000 Hz	1 (Vpp)	2.2 V
1500 Hz	1 (Vpp)	2.2 V

Result: We can see that for different values of frequency in the function generator, we are getting the same value of voltage in the oscilloscope.

B) Oscilloscope connected with DC power supply:


For different values in the function generator, the DC power supply shows the same value (straight line) of voltage in the oscilloscope.



TASK 2:

AC Signal Measurement


- NI compactDAQ 9215 is connected in parallel with the oscilloscope.
- Select values for “sampling rate” and “number of samples” in the LabView for a time varying signal.
- Fixed value for function generator as given,
Frequency = 1000 Hz, Amplitude = 5 V_{pp}




- Below are the graphs which shows for different sampling rate and samples.

GRAPHS:


Rate : 500 Number of Samples : 50



Rate : 3752 Number of Samples : 50




Rate : 3752 Number of Samples : 2000




Rate : 20000

Number of Samples : 2000




Rate : 500


Number of Samples : 2000



Rate : 2000 Number of Samples : 200



Rate : 2000 Number of Samples : 2000




Observation:

- We can see that with increasing sample values, the sine waves formed in the graph increases with respect to time. Whereas when the sample decreases, the wave form decreases.
i.e. **Number of samples** ↑↑(increase), then **Amplitude** ↑↑(increase)
- We also observe that, when the sampling rate increases, the peak to peak value in the graph also increases with respect to time.
i.e. **Sampling rate** ↑↑(increase), then **Peak to peak** ↑↑(increase)


Task 3

Resistance/Power Measurement across a potentiometer w/ Function Generator


- Fixed value for function generator as given,
Frequency = 100 Hz, Amplitude = 3 Vpp
- We can select the value for “sampling rate” and “number of samples” in the LabView.



Function Generator (100Hz
3Vpp)




Rate : 2000 Number of Samples : 2000




Task 4

Find the maximum power output of a function Generator


- The power dissipation is calculated by adding the function in LabView across the two resistors.
- Tune the potentiometer to set the resistance.
- The suggested resistance value of the potentiometer are 20, 40, 60, 80, 100, 150 and 200 ohms.




Potentiometer Resistance = 20 ohm




Potentiometer Resistance = 40 ohm



Potentiometer Resistance = 100 ohm





From the graph, we can see that maximum power occurs at resistance 50 ohms.